

Ultrapeel® II – Microdermabrasion Unit

Easy to use

No clogging

Sterile disposables accessories

Patented 45°, v- shaped, light-weight handpiece

Power adjustment features

Durable



New Features:

Digital crystal flow adjustment

Digital vacuum adjustment

Built-in compressor

Transderm® Compatible

Suitable for:

Facial and body peeling

Aging skin

Dull, congested and dry skin

Enlarged pores

Superficial rhytids

Minor scars

Localized hyperpigmentation

Pre-treatment for further applications (Laser and Chemical Peel)

Mattioli Engineering is proud to introduce our most innovative and technologically advanced medical microdermabrasion system, the Ultrapeel II ®. The Ultrapeel II ® crystal microdermabrasion system is suitable for dermatologists, plastic and surgeons, and medical doctors. The Ultrapeel® II allows users to perform a range of treatments in a safe, hygienic, and sterile environment while avoiding cross contamination risks.

Ultrapeel II ® has a disposable delivery kit that exfoliates the skin through a jet of inert corundum crystals through the patented 45 degree, V-shaped handpiece, while abrading the skin in a gentle but efficient way.

The new crystal flow regulation feature provides users with the opportunity to save crystals. This, in addition to the maximum vacuum power setting, allows the machine to be configured for any type of use and is completely safe for both patients and users.

The Ultrapeel II ® has a built-in compressor for critical applications. The foot pedal that comes with the Ultrapeel II ® enables air compression for deeper treatments by showering the patient's skin with sterile corundum and then vacuuming it out.

The Ultrapeel II ® is the most reliable and versatile machine on the market.

The easy to use accessories, sterilization features, and the power regulations make the Ultrapeel II ® the standard of care in the medical field for any patient and for any skin condition.

The Ultrapeel® was the first microdermabrasion system to be FDA cleared in the United States. FDA cleared - December 1996.

The Ultrapeel II ® is compatible with the Transderm® Ionto transdermal delivery system. This is the only system FDA cleared as the alternative to injections. *(Please review the Transderm Ionto data sheet for details).*

ACCESSORIES AND DISPOSABLES



*Pre-Packed, Pre-Sterilized,
Disposable Delivery Kit,
accessories and corundum
#30000385 Corundum
bottle set
#30010010 Disposable
Tubing kit*



*45 Deg. V-Shaped
Patented disposable
abrading tip
#10010142*

TECHNICAL SPECIFICATIONS

POWER SUPPLY:

Supply Voltage: 110, single phase
Frequency: 60 Hz
Power absorption35VA
Fuses: 2 x 4 AT at 110V, 60Hz
Class (according to IEC 601-1)1 type
BF

PERFORMANCES:

Suction power: from 0.2 to 0.8 bars
(below 10000 mt altitude)

Abrading media:

High purity medical grade corundum
powder

PHYSICAL:

Dimension: 30x30x50 cm
12x12x20 in
Weight: 20 kg (44 lb)

STERILIZATION

Sterilization is of paramount importance when the operator is dealing with different patients and skin conditions.

Cross contamination and debris may affect the patient if the microdermabrader is not perfectly cleaned and sterile.

Only Mattioli Engineering guarantees safe and sterile operations through its disposable accessories and consumables.

Mattioli Engineering - A New Science in Skin Care™

Mattioli Engineering Italia SPA

P.zza di Bellosguardo 11
50124 FIRENZE—ITALY
Phone.: +39 055 221603
Fax: +39 055 221735
info@mattioliengineering.it

**Saratoga Aesthetics**

6909 Sandy Shore Road
Columbia, SC 29206
1.888.785.7546 phone
1.803.788.7776 fax
info@saratogaaesthetics.com
www.SaratogaAesthetics.com