

EUFOTON[®]
medical lasers

endolift[®]

85.5%
Patient
satisfaction*

* 2017-survey based on 2.000 treatments performed.

Endolift[®]

The best non-surgical treatment to boost skin restructuring,
reduce cutaneous laxity and excessive fat.

ONLY ONE TREATMENT | NO ANESTHESIA | IMMEDIATE RESULTS

ENDOLIFT®

ENDOLIFT® is a minimally invasive laser treatment that uses Eufoton® innovative laser **LASEmaR® 1500** (certified and approved by the US FDA for laser assisted liposuction), to stimulate both deep and superficial layers of the skin, tighten and retract the connective septum, stimulate new dermal collagen formation and reduce excessive fat when necessary.

The wavelength of **LASEmaR® 1500** (1470nm) has an ideal interaction with water and fat, which activates neo-collagenesis and the metabolic functions in the extracellular matrix. This results in skin retraction and tightening.

The office-based **ENDOLIFT®** treatment requires Eufoton® specific FTF micro optical fibers, (different calibers depending on the area to treat) that are easily inserted, without any incisions or anesthetic, under the skin directly in the superficial hypodermis, creating a micro-tunnel orientated along the anti-gravitational vectors and, after the treatment, the fibers are removed.

While passing through the dermis, these FTF micro optical fibers act like an intradermal light path and transmit laser energy, offering significant, visible results. The procedure involves minimal to no downtime and it does not have the pain or recovery time that is associated with surgical procedures. Patients can return to work and normal activity within a few hours.

Results are both immediate and long-term. The area will continue to improve for several months following the **ENDOLIFT®** procedure as additional collagen builds in the deep layers of the skin. Results can be enhanced with the combination of **LASEmaR® 1500** and **LIGHTSCAN™** - a fractional, non-ablative laser treatment.

ENDOLIFT® MAIN INDICATIONS

For the areas of initial and middle skin laxity in face and body:

Face

- Lower eyelid
- Mid and lower face
- Mandibular border
- Under chin
- Neck

Body

- Inner arm
- Abdomen and peri-umbilical area
- Inner thigh
- Knee
- Ankle


ENDOLIFT® ADVANTAGES

- Office based procedure
- No anesthesia, just cooling
- Safe and immediate visible results
- Long term effect
- Just one session
- No incisions
- Minimal or no post treatment recovery time
- Possibility of therapeutic combination with fractional non-ablative laser treatment (**LIGHTSCAN™**)

TWO DIFFERENT FTF BARE FIBERS



LOWER EYELID: before - after



CHIN: before - after



MANDIBULAR CONTOUR: before - after



GLUTEUS: before - after



KNEE: before - after



ANKLES: before - after



Before and after photographs are courtesy of D. Cassuto MD, R. Dell'Avanzato MD, R. Forte MD.

EUFOTON®
medical lasers

Eufoton Srl
Via Flavia 23/1 - 34148 Trieste - Italy
Ph. +39 040 9899082 - Fax +39 040 9890555
www.endolift.com - www.eufoton.com
info@eufoton.com

Authorized Distributor in USA:
Saratoga Aesthetics
1.888.785.7546
www.SaratogaAesthetics.com