

“More Effective, Safer and Easier to Use”

Contact Us:

Phone:
(888) 785-7546

Electronic:
info@saratogaesthetics.com




NaturalLight
PULSED LIGHT SYSTEM

**Innovative Aesthetic
Technology**

How NaturaLight Works

Pulsed light allows for light absorption by blood (hemoglobin) and pigment (melanin) for a full spectrum of dermatologic procedures. The NaturaLight will adjust parameters for a specialized treatment for each indication.

Wavelengths

- Ar: 420-1200 nm
- Hair Removal: 640-1200 nm
- Skin Rejuvenation: 530-1200 nm
- NdYag: 1064 nm
- Vascular: 510-1200 nm

Frequency: Up to 1 Hz

Display: LCD Screen

Energy Density: Up to 35 J/cm²

Power Supply: 110 V, 50/60 Hz

Spot Size: 40 X 10 mm

Pulses: 1-8 Pulses

Weight: 50 lbs

Why NaturaLight

- Automatic Parameter Selection
- Variable Pulsing
- Multiple Applications
- Low Cost of Ownership

Types of Treatments

- Hair Removal
- Skin Rejuvenation
- Vascular Lesions
- Nd: YAG Hand Piece



Before & Afters:



BEFORE



AFTER 3 TRMTS



BEFORE



AFTER 4 TRMTS



BEFORE



AFTER 3 TRMTS